



# DIOPSYS® NOVA-ERG

OFFICE BASED PATTERN ELECTRORETINOGRAPHY TESTING



Electrode	Diopsys	DTL	Gold Foil
Position	Skin	Cornea (invasive)	Cornea (invasive)
Repeatability	Good*	Sensitive to eye movements (noise)	Sensitive to eye movements (noise)
Usage	Disposable	Disposable	Disposable
Risk for Corneal Damage and Eye Irritation	Minimal	High	High
Need for Anesthetic	No	Occasionally	Frequently
Application	Easy	Very Difficult	Very Difficult

The easy-to-use Diopsys® NOVA-ERG Lid Sensors are comfortable for the patient and convenient for the doctor.